


AlerTox® Sticks

Rapid tests for allergen detection in food,
drinks, and surfaces for internal Quality Control.



Egg



Crustacean



Hazelnut



β -Lactoglobulin



Fish



Almond



Casein



Mustard Seed



Total Milk



Soy


AlerTox® Sticks advantages:

- » Quick and Reliable
- » Results in 10 minutes
- » No Cross-Reactivity
- » High Sensitivity

AlerTox[®] Sticks are a practical tool for the control of allergen content in raw materials, final products, working surfaces.

Industrial application

- » Semi-quantification of allergens in food
- » Food quality control
- » Trace allergen contamination in food, liquid and surfaces
- » Implementation and fulfillment of HACCP, ISO 22000, BRC, and IFS compliance



Product Name	Ordering No.	Reference No.	Quantity
AlerTox[®] Sticks Egg	KIT3025	KT-5899	10 tests
AlerTox [®] Sticks Crustacean	KIT3036	KT-6126	10 tests
AlerTox[®] Sticks Hazelnut	KIT3035	KT-6129	10 tests
AlerTox [®] Sticks Fish	KIT3038	KT-6131	10 tests
AlerTox[®] Sticks β-Lactoglobulin	KIT3018	KT-5782	10 tests
AlerTox [®] Sticks Almond	KIT3033	KT-6127	10 tests
AlerTox[®] Sticks Mustard Seed	KIT3040	KT-6370	10 tests
AlerTox [®] Sticks Soy	KIT3028	KT-6125	10 tests
AlerTox[®] Sticks Casein	KIT3021	KT-5781	10 tests
AlerTox [®] Sticks Total Milk	KIT3015	KT-6495	10 tests

NOTE: As with any qualitative or semiquantitative test, results should be periodically validated with a quantitative test performed in an outside laboratory.